

Non Sterile Adhesive Laminate Products

Medical Grade Adhesive: Based on Hydrocolloid Adhesive Platform

Features:

- Interacts with skin moisture to provide a secure seal
- Gentle and Secure adhesion that is non irritating to skin
- High fluid absorption - adjustable based on adhesive formula and thickness
- Semi occlusive Adhesive mass plus outer film cover forms barrier to external contaminants
- Silicone release paper - facilitates ease of removal from adhesive and application to skin
- Flexible and conformable - molds to difficult body contours
- No colophany, latex , or animal derived by products
- Product made in USA
- Manufacturing Plant: FDA - approved facility
ISO 13485 certified

Uses:

- Professional Wound dressings
- Ostomy adhesive barriers
- Incontinence Fixation adhesives
- Consumer (OTC) first aid dressings
- Consumer heel blister products
- Consumer blemish covers

Standard Roll/Sheet Data

Dimensions:	Roll Goods:	Standard Width:	240MM
		Standard Length	100 meters
		Thickness: Varied :	i.e. 0.2 - 1MM
	Sheets	Standard Width:	240 MM
		Thickness: Varied :	i.e. 0.2 - 3MM

- Release paper :** silicone with variable number of release liner slits
- Film Covers:** Various : Polyethylene, Polyurethane, non woven, or custom
- Adhesives:** Patented SureSkin Hydrocolloid or custom variations
- Packaging:** Rolls are wrapped in plastic and boxed
- Artwork:** Print on release liner and some cover films can be customized

EuroMed Laminate Products

Better Performance...Better Value....Better by design




411 Clinton Avenue
Northvale, NJ 07647
Contact:
Steven Hartnack Vice President of Sales
shartnack@euromedinc.com
P. (201) 750-1840 x 11
USA only 1-877 EURO 595
F. (201) 750-1841
www.euromedinc.com





Better by Design ●

Laminate Adhesive Roll Good Product Data:

Relative values based on average roll thickness of 0.5 mm



-  H014 Roll Good Non-Sterile
-  H020 Roll Good Non-Sterile

